



# SUPER VIEW® PIVOT

for use with the SUPER VIEW® and BIOM® wide angle viewing systems




## See what you have been missing with the disposable PIVOT™ Lens Set for vitreo-retinal surgery

130° HD Wide Field View



72° Ultra HD View



SUPER VIEW® PIVOT

- Effortlessly rotate between 130° HD Wide Field peripheral and 72° Ultra HD 3D macular views
- Convenience of two lenses in one disposable set
- Single use for a crisp clear view every time
- Increased working distance for instrument manipulation and to minimize fogging
- Outstanding depth of field
- No assembly required
- Sets available with the HTC™ (Hassan-Tornambe Contact) Lens



**Please contact the customer service department for pricing**

Insight Instruments, Inc. | 2580 S.E. Willoughby Blvd. Stuart, FL 34994  
Tel: 800-255-8354 (US Only) | 772-219-9393 | Fax: 772-219-9342  
[www.insightinstruments.com](http://www.insightinstruments.com) | [CustomerService@insightinstruments.com](mailto:CustomerService@insightinstruments.com)